



Product Catalogue 产品目录

网址: www.scienova.com

● about scienova

关于 scienova

scienova GmbH is located at the research and innovation center Jena in Thuringia (Germany).

Since its foundation in the year 2008, scienova GmbH has been focusing on the development and distribution of automated dialysis cartridges, which make sample dialysis for research and diagnostic processes more productive.

scienova GmbH was founded by scientists from the Institute of Biochemistry at the University Hospital Jena and has built up a large network of partners and distributors since then.

Innovative companies and R&D facilities in the fields of biotech technology, medicine, pharmacy and diagnostics are among our customers.

Due to our small organizational structure, we are able to respond flexibly and very quickly to customer-specific inquiries and can implement individual products with parameters according to your specifications.

We invite you to get to know us and our products and to challenge us with your problems.

scienova 有限公司坐落在德意志联邦共和国图林根州耶拿市研究创业基地。

自从2008年公司成立以来，scienova公司一直致力于研发和经营自动化透析滤芯产品，大大提高了研发和诊断过程所需的样品透析效率。


scienova公司由Jena大学附属医院生物化学研究所科学家团队创建，现已建立具有完善的合作和经营网络，并为生物技术、医药和医学检验行业领域的创业公司和研发机构提供服务。

依靠一种小型公司具备的组织和结构优势，我们能够灵活、快捷地应对不断变化的客户需求，按照客户提出的技术参数，提供个体化产品服务。

我们诚邀您，认识我们，了解我们的产品，向我们提出问题，挑战我们。

Xpress Dialyzer Family / Xpress 透析器家

	MD100	ED300	MD1000	MD300	XDR-1
Sample volume 样品体积	10-100 µl	50-300 µl	150-1000 µl	50-300 µl	1-9 ml
Deep-well plates 深孔板	96-well 96-孔板	96-well 96-孔板	48-well 48-孔板	96-well 96-孔板	6-well 6-孔板
Membrane material 膜材料	Regenerated cellulose 再生纤维素				
Membrane pore sizes 膜孔径大小	2 3.5 6-8 12-14 20 kDa				
Temperature 温度	1-60 °C				
pH value pH值	4-8				
Membrane material 膜材料	Bio-Xell® 再生纤维素				
Membrane pore sizes 膜孔径大小	140 kDa				
Temperature 温度	1-60 °C				
pH value pH值	5-8				



处理视频




goo.gl/jnROnh

Xpress Dialyzer / Xpress 透析器

MD100

Xpress Micro Dialyzer / Xpress 微透析器




Article 产品	MWCO 截留分子量	Number of samples / Description 样品数量/说明	Article number 产品型号
MD100 (10-100 µl)	3.5 kDa	96 pcs. in deep well plate (12 cartridges), sealing sheet 96孔一板（12条滤芯），密封盖	40072
		8 pcs. (1 cartridge) 8 孔（1条滤芯）	40071
		12 pcs. (segments) 12孔（分格）	40782
		48 pcs. (6 cartridges), GridKit48, deep well plate incl. grid and lid 48孔（6条滤芯），48孔板，含网格和密封盖	40840
		60 pcs. (segments, incl. universal rack, beaker, stirrer and forceps, float) 60孔（格子，含多功能架子，拆包器，搅拌棒，镊子和平板玻璃）	40071-X60
		288 pcs. (segments, incl. universal rack, beaker, stirrer and forceps, float) 288孔（格子，含多功能架子，拆包器，搅拌棒，镊子和平板玻璃）	40071-X288
	6-8 kDa	96 pcs. in deep well plate (12 cartridges), sealing sheet 96孔板（12条滤芯），密封盖	40076
		8 pcs. (1 cartridge) 8 孔（1条滤芯）	40075
		12 pcs. (segments) 12孔（分格）	40783
		48 pcs. (6 cartridges), GridKit48, deep well plate incl. grid and lid 48孔（6条滤芯），48孔板，含网格和密封盖	40841
		60 pcs. (segments, incl. universal rack, beaker, stirrer and forceps) 60孔（格子，含多功能架子，拆包器，搅拌棒，镊子和平板玻璃）	40075-X60
		288 pcs. (segments, incl. universal rack beaker, stirrer and forceps) 288孔（格子，含多功能架子，拆包器，搅拌棒，镊子，平板玻璃）	40075-X288
	12-14 kDa	96 pcs. in deep well plate (12 cartridges), sealing sheet 96孔板（12条滤芯），密封盖	40078
		8 pcs. (1 cartridge) 8 孔（1条滤芯）	40077
		12 pcs. (segments) 12孔（分格）	40784
		48 pcs. (6 cartridges), GridKit48, deep well plate incl. grid and lid 48孔（6条滤芯），48孔板，含网格和密封盖	40842
		60 pcs. (segments, incl. universal rack, beaker, stirrer and forceps, float) 60孔（格子，含多功能架子，拆包器，搅拌棒，镊子和平板玻璃）	40077-X60
		288 pcs. (segments, incl. universal rack, beaker, stirrer and forceps, float) 288孔（格子，含多功能架子，拆包器，搅拌棒，镊子，平板玻璃）	40077-X288

Xpress Dialyzer / Xpress 透析器

Article 产品	MWCO 截留分子量	Number of samples / Description 样品数量/说明	Article number 产品型号
MD100 (10-100 µl)	2 kDa	96 pcs. in deep well plate (12 cartridges), sealing sheet 96孔板 (12条滤芯), 密封盖	40932
		8 pcs. (1 cartridge) 8 孔 (1条滤芯)	40933
	20 kDa	96 pcs. in deep well plate (12 cartridges), sealing sheet 96孔板 (12条滤芯), 密封盖	40934
		8 pcs. (1 cartridge) 8 孔 (1条滤芯)	40935
	140 kDa	96 pcs. in deep well plate (12 cartridges), sealing sheet 96孔板 (12条滤芯), 密封盖	40930
	Bio-Xell® (biocel- lulose)	8 pcs. (1 cartridge) 8 孔 (1条滤芯)	40931

MD300

Xpress Micro Dialyzer / Xpress 微透析器



Article 产品	MWCO 截留分子量	Number of samples / Description 样品数量/说明	Article number 产品型号
MD300 (50-300 µl)	3.5 kDa	96 pcs. in deep well plate (12 cartridges), sealing sheet 96孔板 (12条滤芯), 密封盖	40787
		8 pcs. (1 cartridge) 8 孔 (1条滤芯)	40786
		48 pcs. (6 cartridges), GridKit48, deep well plate incl. grid and lid 48孔 (6条滤芯), 48孔板, 含网格和密封盖	40843
		12 pcs. (segments) 12孔 (分格)	40786-X12
	6-8 kDa	96 pcs. in deep well plate (12 cartridges), sealing sheet 96孔板 (12条滤芯), 密封盖	40789
		8 pcs. (1 cartridge) 8 孔 (1条滤芯)	40788
		48 pcs. (6 cartridges), GridKit48, deep well plate incl. grid and lid 48孔 (6条滤芯), 48孔板, 含网格和密封盖	40844
		12 pcs. (segments) 12孔 (分格)	40788-X12

Xpress Dialyzer / Xpress 透析器

Article 产品	MWCO 截留分子量	Number of samples / Description 样品数量/说明	Article number 产品型号
MD300 (50-300 µl)	12-14 kDa	96 pcs. in deep well plate (12 cartridges), sealing sheet 96孔板 (12条滤芯), 密封盖	40791
		8 pcs. (1 cartridge) 8 孔 (1条滤芯)	40790
		48 pcs. (6 cartridges), GridKit48, deep well plate incl. grid and lid 48孔 (6条滤芯), 48孔板, 含网格和密封盖	40845
		12 pcs. (segments) 12孔 (分格)	40790-X12
	2 kDa	96 pcs. in deep well plate (12 cartridges), sealing sheet 96孔板 (12条滤芯), 密封盖	40961
	20 kDa	8 pcs. (1 cartridge) 8 孔 (1条滤芯)	40962
		96 pcs. in deep well plate (12 cartridges), sealing sheet 96孔板 (12条滤芯), 密封盖	40963
	140 kDa	8 pcs. (1 cartridge) 8 孔 (1条滤芯)	40964
		96 pcs. in deep well plate (12 cartridges), sealing sheet 96孔板 (12条滤芯), 密封盖	40965
Bio-Xell® (biocel- lulose)		8 pcs. (1 cartridge) 8 孔 (1条滤芯)	40966

GridKit48
Liquid Handling
Position Grid
液体处理位置网格



Xpress Dialyzer / Xpress 透析器

ED300

Xpress Equilibrium Dialyzer / Xpress 平衡透析器



Article 产品	MWCO 截留分子量	Number of samples / Description 样品数量/说明	Article number 产品型号
ED300 (50 - 300 µl)	3.5 kDa	96 pcs. in deep well plate (12 cartridges), sealing sheet 96孔板 (12条滤芯), 密封盖	40891
		8 pcs. (1 cartridge) 8 孔 (1条滤芯)	40890
	6-8 kDa	96 pcs. in deep well plate (12 cartridges), sealing sheet 96孔板 (12条滤芯), 密封盖	40893
		8 pcs. (1 cartridge) 8 孔 (1条滤芯)	40892
	12-14 kDa	96 pcs. in deep well plate (12 cartridges), sealing sheet 96孔板 (12条滤芯), 密封盖	40895
		8 pcs. (1 cartridge) 8 孔 (1条滤芯)	40894
	2 kDa	96 pcs. in deep well plate (12 cartridges), sealing sheet 96孔板 (12条滤芯), 密封盖	40957
		8 pcs. (1 cartridge) 8 孔 (1条滤芯)	40956
	20 kDa	96 pcs. in deep well plate (12 cartridges), sealing sheet 96孔板 (12条滤芯), 密封盖	40939
		8 pcs. (1 cartridge) 8 孔 (1条滤芯)	40938
	140 kDa	96 pcs. in deep well plate (12 cartridges), sealing sheet 96孔板 (12条滤芯), 密封盖	40937
		8 pcs. (1 cartridge) 8 孔 (1条滤芯)	40936
Bio-Xell® (biocellulose)		8 pcs. (1 cartridge) 8 孔 (1条滤芯)	40936

Two-channel system:

Pipette in sample and remove test dialysate from the top of the Dialyzer without removing sample cartridges.

双通道系统:

在未拆卸滤芯（膜）条件下，可以从透析器头部注入样品或吸取透析液。



Xpress Dialyzer / Xpress 透析器

MD1000

Xpress Mini Dialyzer / Xpress 小型透析器



Article 产品	MWCO 截留分子量	Number of samples / Description 样品数量/说明	Article number 产品型号
MD1000 (50-1000 µl)	3.5 kDa	6 segments in tubes, float 离心管式浮动滤芯格， 6只	40870
	6-8 kDa	6 segments in tubes, float 离心管式浮动滤芯格， 6只	40873
	12-14 kDa	6 segments in tubes, float 离心管式浮动滤芯格， 6只	40876
	2 kDa	6 segments in tubes, float 离心管式浮动滤芯格， 6只	40991
	20 kDa	6 segments in tubes, float 离心管式浮动滤芯格， 6只	40994
	140 kDa Bio-Xell®	6 segments in tubes, float 离心管式浮动滤芯格， 6只	40997
	3.5 kDa	6 segments in tubes (plus 6 additional tubes), float 离心管式浮动滤芯格（附加6只空心离心管）， 6只	40871
	6-8 kDa	6 segments in tubes (plus 6 additional tubes), float 离心管式浮动滤芯格（附加6只空心离心管）， 6只	40874
	12-14 kDa	6 segments in tubes (plus 6 additional tubes), float 离心管式浮动滤芯格（附加6只空心离心管）， 6只	40877
	2 kDa	6 segments in tubes (plus 6 additional tubes), float 离心管式浮动滤芯格（附加6只空心离心管）， 6只	40992
	20 kDa	6 segments in tubes (plus 6 additional tubes), float 离心管式浮动滤芯格（附加6只空心离心管）， 6只	40995
	140 kDa Bio-Xell®	6 segments in tubes (plus 6 additional tubes), float 离心管式浮动滤芯格（附加6只空心离心管）， 6只	40998
	3.5 kDa	48 segments in 48-deep well plate (Kit), sealing sheet 48孔板式滤芯格（盒子）， 密封盖	40872
	6-8 kDa	48 segments in 48-deep well plate (Kit), sealing sheet 48孔板式滤芯格（盒子）， 密封盖	40875
	12-14 kDa	48 segments in 48-deep well plate (Kit), sealing sheet 48孔板式滤芯格（盒子）， 密封盖	40878
	2 kDa	48 segments in 48-deep well plate (Kit), sealing sheet 48孔板式滤芯格（盒子）， 密封盖	40993
	20 kDa	48 segments in 48-deep well plate (Kit), sealing sheet 48孔板式滤芯格（盒子）， 密封盖	40996
	140 kDa Bio-Xell®	48 segments in 48-deep well plate (Kit), sealing sheet 48孔板式滤芯格（盒子）， 密封盖	40999

Xpress Dialyzer / Xpress 透析器


XDR-1

Xpress Dialysis Reactor / Xpress 透析反应器



Article 产品	MWCO 截留分子量	Number of samples / Description 样品数量/说明	Article number 产品型号
XDR-1 (1-19 ml)	3.5 kDa	1 x Dialyzer, 1 x floater, 3 x plug 1只透析器, 1只浮动框, 3 塞栓	41023
	6-8 kDa	1 x Dialyzer, 1 x floater, 3 x plug 1只透析器, 1只浮动框, 3 塞栓	41026
	12-14 kDa	1 x Dialyzer, 1 x floater, 3 x plug 1只透析器, 1只浮动框, 3 塞栓	41029
	3.5 kDa	6 x Dialyzer, 18 x plug, 6 x disposable syringe, sealing sheet, 6-deep well plate (80 ml) 6只透析器, 18只塞栓, 6只一次性注射器, 密封盖, 6-孔板 (80 ml)	41024
	6-8 kDa	6 x Dialyzer, 18 x plug, 6 x disposable syringe, sealing sheet, 6-deep well plate (80 ml) 6只透析器, 18只塞栓, 6只一次性注射器, 密封盖, 6-孔板 (80 ml)	41027
	12-14 kDa	6 x Dialyzer, 18 x plug, 6 x disposable syringe, sealing sheet, 6-deep well plate (80 ml) 6只透析器, 18只塞栓, 6只一次性注射器, 密封盖, 6-孔板 (80 ml)	41030
	3.5 kDa	6 x Dialyzer, 3 x floater, 18 x plug, 6 x disposable syringe 6只透析器, 3只浮动框, 18只塞栓, 6只一次性注射器	41025
	6-8 kDa	6 x Dialyzer, 3 x floater, 18 x plug, 6 x disposable syringe 6只透析器, 3只浮动框, 18只塞栓, 6只一次性注射器	41028
	12-14 kDa	6 x Dialyzer, 3 x floater, 18 x plug, 6 x disposable syringe 6只透析器, 3只浮动框, 18只塞栓, 6只一次性注射器	41031

Article 产品	Number of samples / Description 样品数量/说明	Article number 产品型号
XDR Accessories	2 x tube connection (angle connection) 离心管联结 (角型联结), 2条	41038
	2 x plug small 2只塞栓, 小	41039
	1 x plug big 1只塞栓, 大	41040
	2 x straight connection 2条直接联结	41041
	1 x silicone tube (1.5 m) 硅胶管 (1.5 m), 1只	41042
	1 x 6-deep well plate, 1 x sealing sheet 1只6-孔板, 1只密封盖	41046
	6 x disposable syringe (20 ml) 6只一次性注射器 (20 ml)	41044



GridKit48

Liquid Handling Position Grid

48-孔板液体工作平台

- Available with scienova Xpress Dialyzer MD100, MD300 and ED300
- Free positions in grid can be used for buffer handling while dialysis
- Designed and tested as automated liquid handling device



- 此为scienova Xpress 透析器MD100, MD300和ED300试剂盒提供的附加产品
- 透析期间, 该工作平台的空位可以用于缓冲液操作
- 为液体操作设计和制作的自动化工作平台

● our quality 我们的质量保障

scienova GmbH has been certified for the quality management system according to ISO 9001:2015 since the year 2016.

We guarantee quality assurance throughout the entire manufacturing process.

We ensure the consistent quality of our products through regular quality controls and document these quality data in batch certificates.

With our longtime market experience and scientific know-how we pursue a continuous product development process. Our focus is set on the improvement in dialysis efficiency by innovative high-performance products.

2016年以来, scienova 生物技术拥有ISO9001:2015 质量管理体系资质。我们为整体生产过程提供质量安全保障。

我们产品质量的稳定性体现在定期质量控制和批次质检合格证。

我们拥有长期积累的市场经验和科技思路, 能够持续不断地开展产品研发活动。我们利用新兴的高效科技产品, 不断更新透析效率是我们的核心竞争力所在。

Xpress Dialyzer / Xpress 透析器

Mixing Box

Xpress Magnetic Mixing Dialysis Box / Xpress 快速磁铁混匀透析槽



(magnetic stirrer is not part of the product)

(本产品不含磁铁胶棒)

Article 产品	MWCO 截留分子量	Number of samples / Description 样品数量/说明	Article number 产品型号
Xpress Magnetic Mixing Dialysis Box MD100 (10-100 µl)	3.5 kDa	96 samples (12 cartridges MD100), dialysis and handling box, magnetic stirrer, sealing film, self adhesive positioning bumpons 96份样品（12条滤芯MD100），透析和作业槽，磁铁胶棒，密封膜，粘扣带	40940
		Refill kit, 96 samples (12 cartridges MD100), handling box, sealing film 填充试剂盒，96份样品（12条滤芯MD100），作业槽，密封膜	40941
	6-8 kDa	96 samples (12 cartridges MD100), dialysis and handling box, magnetic stirrer, sealing film, self adhesive positioning bumpons 96份样品（12条滤芯MD100），透析和作业槽，磁铁胶棒，密封膜，粘扣带	40942
		Refill kit, 96 samples (12 cartridges MD100), handling box, sealing film 填充试剂盒，96份样品（12条滤芯MD100），作业槽，密封膜	40943
	12-14 kDa	96 samples (12 cartridges MD100), dialysis and handling box, magnetic stirrer, sealing film, self adhesive positioning bumpons 96份样品（12条滤芯MD100），透析和作业槽，磁铁胶棒，密封膜，粘扣带	40944
		Refill kit, 96 samples (12 cartridges MD100), handling box, sealing film 填充试剂盒，96份样品（12条滤芯MD100），作业槽，密封膜	40945
	2 kDa	96 samples (12 cartridges MD100), dialysis and handling box, magnetic stirrer, sealing film, self adhesive positioning bumpons 96份样品（12条滤芯MD100），透析和作业槽，磁铁胶棒，密封膜，粘扣带	40946
		Refill kit, 96 samples (12 cartridges MD100), handling box, sealing film 填充试剂盒，96份样品（12条滤芯MD100），作业槽，密封膜	40947
	20 kDa	96 samples (12 cartridges MD100), dialysis and handling box, magnetic stirrer, sealing film, self adhesive positioning bumpons 96份样品（12条滤芯MD100），透析和作业槽，磁铁胶棒，密封膜，粘扣带	40948
		Refill kit, 96 samples (12 cartridges MD100), handling box, sealing film 填充试剂盒，96份样品（12条滤芯MD100），作业槽，密封膜	40958
	140 kDa	96 samples (12 cartridges MD100), dialysis and handling box, magnetic stirrer, sealing film, self adhesive positioning bumpons 96份样品（12条滤芯MD100），透析和作业槽，磁铁胶棒，密封膜，粘扣带	40959
	Bio-Xell®	Refill kit, 96 samples (12 cartridges MD100), handling box, sealing film 填充试剂盒，96份样品（12条滤芯MD100），作业槽，密封膜	40960

Xpress Dialyzer / Xpress 透析器

Article 产品	MWCO 截留分子量	Number of samples / Description 样品数量/说明	Article number 产品型号
Xpress Magnetic Mixing Dialysis Box MD300 (50-300 µl)	3.5 kDa	96 samples (12 cartridges MD300), dialysis and handling box, magnetic stirrer, sealing film, self adhesive positioning bumpons 96份样品（12条滤芯MD300），透析和作业槽，磁铁胶棒，密封膜，粘扣带	40969
		Refill kit, 96 samples (12 cartridges MD300), handling box, sealing film 填充试剂盒，96份样品（12条滤芯MD300），作业槽，密封膜	40970
	6-8 kDa	96 samples (12 cartridges MD300), dialysis and handling box, magnetic stirrer, sealing film, self adhesive positioning bumpons 96份样品（12条滤芯MD300），透析和作业槽，磁铁胶棒，密封膜，粘扣带	40971
		Refill kit, 96 samples (12 cartridges MD300), handling box, sealing film 填充试剂盒，96份样品（12条滤芯MD300），作业槽，密封膜	40972
	12-14 kDa	96 samples (12 cartridges MD300), dialysis and handling box, magnetic stirrer, sealing film, self adhesive positioning bumpons 96份样品（12条滤芯MD300），透析和作业槽，磁铁胶棒，密封膜，粘扣带	40973
		Refill kit, 96 samples (12 cartridges MD300), handling box, sealing film 填充试剂盒，96份样品（12条滤芯MD300），作业槽，密封膜	40974
	2 kDa	96 samples (12 cartridges MD300), dialysis and handling box, magnetic stirrer, sealing film, self adhesive positioning bumpons 96份样品（12条滤芯MD300），透析和作业槽，磁铁胶棒，密封膜，粘扣带	41003
		Refill kit, 96 samples (12 cartridges MD300), handling box, sealing film 填充试剂盒，96份样品（12条滤芯MD300），作业槽，密封膜	41004
	20 kDa	96 samples (12 cartridges MD300), dialysis and handling box, magnetic stirrer, sealing film, self adhesive positioning bumpons 96份样品（12条滤芯MD300），透析和作业槽，磁铁胶棒，密封膜，粘扣带	40975
		Refill kit, 96 samples (12 cartridges MD300), handling box, sealing film 填充试剂盒，96份样品（12条滤芯MD300），作业槽，密封膜	40976
	140 kDa	96 samples (12 cartridges MD300), dialysis and handling box, magnetic stirrer, sealing film, self adhesive positioning bumpons 96份样品（12条滤芯MD300），透析和作业槽，磁铁胶棒，密封膜，粘扣带	40977
	Bio-Xell®	Refill kit, 96 samples (12 cartridges MD300), handling box, sealing film 填充试剂盒，96份样品（12条滤芯MD300），作业槽，密封膜	40978
Xpress Magnetic Mixing Dialysis Box ED300 (50-300 µl)	3.5 kDa	96 samples (12 cartridges ED300), dialysis and handling box, magnetic stirrer, sealing film, self adhesive positioning bumpons 96份样品（12条滤芯ED300），透析和作业槽，磁铁胶棒，密封膜，粘扣带	40981
		Refill kit, 96 samples (12 cartridges ED300), handling box, sealing film 填充试剂盒，96份样品（12条滤芯ED300），作业槽，密封膜	40982
	6-8 kDa	96 samples (12 cartridges ED300), dialysis and handling box, magnetic stirrer, sealing film, self adhesive positioning bumpons 96份样品（12条滤芯ED300），透析和作业槽，磁铁胶棒，密封膜，粘扣带	40983
		Refill kit, 96 samples (12 cartridges ED300), handling box, sealing film 填充试剂盒，96份样品（12条滤芯ED300），作业槽，密封膜	40984

Xpress Dialyzer / Xpress 透析器

Article 产品	MWCO 截留分子量	Number of samples / Description 样品数量/说明	Article number 产品型号
	12-14 kDa	96 samples (12 cartridges ED300), dialysis and handling box, magnetic stirrer, sealing film, self adhesive positioning bumpons 96份样品（12条滤芯ED300），透析和作业槽，磁铁胶棒，密封膜，粘扣带 Refill kit, 96 samples (12 cartridges ED300), handling box, sealing film 填充试剂盒，96份样品（12条滤芯ED300），作业槽，密封膜	40985 40986
	2 kDa	96 samples (12 cartridges ED300), dialysis and handling box, magnetic stirrer, sealing film, self adhesive positioning bumpons 96份样品（12条滤芯ED300），透析和作业槽，磁铁胶棒，密封膜，粘扣带 Refill kit, 96 samples (12 cartridges ED300), handling box, sealing film 填充试剂盒，96份样品（12条滤芯ED300），作业槽，密封膜	40979 40980
	20 kDa	96 samples (12 cartridges ED300), dialysis and handling box, magnetic stirrer, sealing film, self adhesive positioning bumpons 96份样品（12条滤芯ED300），透析和作业槽，磁铁胶棒，密封膜，粘扣带 Refill kit, 96 samples (12 cartridges ED300), handling box, sealing film 填充试剂盒，96份样品（12条滤芯ED300），作业槽，密封膜	40987 40988
	140 kDa	96 samples (12 cartridges ED300), dialysis and handling box, magnetic stirrer, sealing film, self adhesive positioning bumpons 96份样品（12条滤芯ED300），透析和作业槽，磁铁胶棒，密封膜，粘扣带	40989
	Bio-Xell®	Refill kit, 96 samples (12 cartridges ED300), handling box, sealing film 填充试剂盒，96份样品（12条滤芯ED300），作业槽，密封膜	40990

处理视频



goo.gl/jnROnh



Dialysis Tubes / 透析管



Article 产品	MWCO 截留分子量	Samples 样品数量	Description 说明	Article number 产品型号
Mini / 小型 (10-250 µl)	6-8 kDa	10	GeBAflex-tube kit (10) Mini, 6-8 kDa, Float, Manual 小型GeBAflex透析管 (10), 6-8 kDa, 浮动式, 附有说明书	40176
		30	GeBAflex-tube kit (30) Mini, 6-8 kDa, Float, Manual 小型GeBAflex透析管 (30), 6-8 kDa, 浮动式, 附有说明书	40178
	12-14 kDa	10	GeBAflex-tube kit (10) Mini, 12-14 kDa, Float, Manual 小型GeBAflex透析管 (10), 12-14 kDa, 浮动式, 附有说明书	40177
		30	GeBAflex-tube kit (30) Mini, 12-14 kDa, Float, Manual 小型GeBAflex透析管 (30), 6-8 kDa, 浮动式, 附有说明书	40179
Midi / 中型 (100-800 µl)	25 kDa	5	GeBAflex-tube kit (5) Mini, 25 kDa, Float, Manual 小型GeBAflex透析管 (5), 25 kDa, 浮动式, 附有说明书	40899
	1 kDa	5	GeBAflex-tube kit (5) Midi, 1 kDa, Float, Manual 小型GeBAflex透析管 (5), 1 kDa, 浮动式, 附有说明书	40898
	3.5 kDa	10	GeBAflex-tube kit (10) Midi, 3.5 kDa, Float, Manual 中型GeBAflex透析管 (10), 3.5 kDa, 浮动式, 附有说明书	40180
		30	GeBAflex-tube kit (30) Midi, 3.5 kDa, Float, Manual 中型GeBAflex透析管 (30), 3.5 kDa, 浮动式, 附有说明书	40181
	6-8 kDa	10	GeBAflex-tube kit (10) Midi, 6-8 kDa, Float, Manual 中型GeBAflex透析管 (10), 6-8 kDa, 浮动式, 附有说明书	40182
		30	GeBAflex-tube kit (30) Midi, 6-8 kDa, Float, Manual 中型GeBAflex透析管 (30), 6-8 kDa, 浮动式, 附有说明书	40183

Dialysis Tubes / 透析管

Article 产品	MWCO 截留分子量	Samples 样品数量	Description 说明	Article number 产品型号
Maxi / 大型 (0.1-3 ml)	3.5 kDa	5	GeBAflex-tube kit (5) Maxi, 3.5 kDa, Float, Manual 大型GeBAflex透析管 (5), 3.5 kDa, 浮动式, 附有说明书	40184
		15	GeBAflex-tube kit (15) Maxi, 3.5 kDa, Float, Manual 大型GeBAflex透析管 (15), 3.5 kDa, 浮动式, 附有说明书	40185
	6-8 kDa	5	GeBAflex-tube kit (5) Maxi, 6-8 kDa, Float, Manual 大型GeBAflex透析管 (5), 6-8 kDa, 浮动式, 附有说明书	40186
		15	GeBAflex-tube kit (15) Maxi, 6-8 kDa, Float, Manual 大型GeBAflex透析管 (15), 6-8 kDa, 浮动式, 附有说明书	40187
	12-14 kDa	5	GeBAflex-tube kit (5) Maxi, 12-14 kDa, Float, Manual 大型GeBAflex透析管 (5), 12-14 kDa, 浮动式, 附有说明书	40188
		15	GeBAflex-tube kit (15) Maxi, 12-14 kDa, Float, Manual 大型GeBAflex透析管 (15), 12-14 kDa, 浮动式, 附有说明书	40189
	25 kDa	5	GeBAflex-tube kit (5) Maxi, 25 kDa, Float, Manual 大型GeBAflex透析管 (5), 25 kDa, 浮动式, 附有说明书	40896
	50 kDa	15	GeBAflex-tube kit (15) Maxi, 50 kDa, Float, Manual 大型GeBAflex透析管 (15), 50 kDa, 浮动式, 附有说明书	40897
Mega / 超大型 (10-20 ml)	1 KDa	5	GeBAflex-tube kit (5) Mega, 1 kDa, Float, Manual, 10 ml 超大型GeBAflex透析管 (5), 1 kDa, 浮动式, 附有说明书, 10 ml	41000
		5	GeBAflex-tube kit (5) Mega, 1 kDa, Float, Manual, 15 ml 超大型GeBAflex透析管 (5), 1 kDa, 浮动式, 附有说明书, 15 ml	41001
		5	GeBAflex-tube kit (5) Mega, 1 kDa, Float, Manual, 20 ml 超大型GeBAflex透析管 (5), 1 kDa, 浮动式, 附有说明书, 20 ml	41002
	3.5 kDa	10	GeBAflex-tube kit (10) Mega, 3.5 kDa, Float, Manual, 10 ml 超大型GeBAflex透析管 (10), 3.5 kDa, 浮动式, 附有说明书, 10 ml	40190
		10	GeBAflex-tube kit (10) Mega, 3.5 kDa, Float, Manual, 15 ml 超大型GeBAflex透析管 (10), 3.5 kDa, 浮动式, 附有说明书, 15 ml	40193
		10	GeBAflex-tube kit (10) Mega, 3.5 kDa, Float, Manual, 20 ml 超大型GeBAflex透析管 (10), 3.5 kDa, 浮动式, 附有说明书, 20 ml	40196



Dialysis Tubes / 透析管

Article 产品	MWCO 截留分子量	Samples 样品数量	Description 说明	Article number 产品型号
	6-8 kDa	10	GeBAflex-tube kit (10) Mega, 6-8 kDa, Float, Manual, 10 ml 超大型GeBAflex透析管 (10), 6-8 kDa, 浮动式, 附有说明书, 10 ml	40191
		10	GeBAflex-tube kit (10) Mega, 6-8 kDa, Float, Manual, 15 ml 超大型GeBAflex透析管 (10), 6-8 kDa, 浮动式, 附有说明书, 15 ml	40194
		10	GeBAflex-tube kit (10) Mega, 6-8 kDa, Float, Manual, 20 ml 超大型GeBAflex透析管 (10), 6-8 kDa, 浮动式, 附有说明书, 20 ml	40197
	12-14 kDa	10	GeBAflex-tube kit (10) Mega, 12-14 kDa, Float, Manual, 10 ml 超大型GeBAflex透析管 (10), 12-14 kDa, 浮动式, 附有说明书, 10 ml	40192
		10	GeBAflex-tube kit (10) Mega, 12-14 kDa, Float, Manual, 15 ml 超大型GeBAflex透析管 (10), 12-14 kDa, 浮动式, 附有说明书, 15 ml	40195
		10	GeBAflex-tube kit (10) Mega, 12-14 kDa, Float, Manual, 20 ml 超大型GeBAflex透析管 (10), 12-14 kDa, 浮动式, 附有说明书, 20 ml	40198

ZelluTrans



Article 产品	MWCO 截留分子量	width l Ø [mm] 宽 Ø [mm]	ml/ cm	Description 说明	Article number 产品型号
T1	3.5 kDa	19 l 12.1	1.2	Dialysis membrane 5 m, Manual, data sheet 5 m 透析膜, 说明书, 数据表	40202-5
				Dialysis membrane 30 m, Manual, data sheet, scissors and two clamps 30 m 透析膜, 说明书, 数据表, 剪刀, 2只夹子	40202
		46 l 29.3	6.7	Dialysis membrane 5 m, Manual, data sheet 5 m 透析膜, 说明书, 数据表	40203-5
				Dialysis membrane 30 m, Manual, data sheet, scissors and two clamps 30 m 透析膜, 说明书, 数据表, 剪刀, 2只夹子	40203

Dialysis Tubes

Article 产品	MWCO 截留分子量	width l Ø [mm] 宽 Ø [mm]	ml/ cm	Description 说明	Article number 产品型号
T2	6-8 kDa	25 l 16	2.0	Dialysis membrane 5 m, Manual, data sheet 5 m 透析膜, 说明书, 数据表	40213-5
				Dialysis membrane 30 m, Manual, data sheet, scissors and two clamps 30 m 透析膜, 说明书, 数据表, 剪刀, 2只夹子	40213
		32 l 20.4	3.3	Dialysis membrane 5 m, Manual, data sheet 5 m 透析膜, 说明书, 数据表	40214-5
				Dialysis membrane 30 m, Manual, data sheet, scissors and two clamps 30 m 透析膜, 说明书, 数据表, 剪刀, 2只夹子	40214
		40 l 25.5	5.1	Dialysis membrane 5 m, Manual, data sheet 5 m 透析膜, 说明书, 数据表	40215-5
				Dialysis membrane 30 m, Manual, data sheet, scissors and two clamps 30 m 透析膜, 说明书, 数据表, 剪刀, 2只夹子	40215
		50 l 31.9	8.0	Dialysis membrane 5 m, Manual, data sheet 5 m 透析膜, 说明书, 数据表	40216-5
				Dialysis membrane 30 m, Manual, data sheet, scissors and two clamps 30 m 透析膜, 说明书, 数据表, 剪刀, 2只夹子	40216
		25 l 15.9	2.0	Dialysis membrane 5 m, Manual, data sheet 5 m 透析膜, 说明书, 数据表	40222-5
				Dialysis membrane 30 m, Manual, data sheet, scissors and two clamps 30 m 透析膜, 说明书, 数据表, 剪刀, 2只夹子	40222
		33 l 21	3.5	Dialysis membrane 5 m, Manual, data sheet 5 m 透析膜, 说明书, 数据表	40223-5
				Dialysis membrane 30 m, Manual, data sheet, scissors and two clamps 30 m 透析膜, 说明书, 数据表, 剪刀, 2只夹子	40223
T3	12-14 kDa	45 l 28.6	6.4	Dialysis membrane 5 m, Manual, data sheet 5 m 透析膜, 说明书, 数据表	40224-5
				Dialysis membrane 30 m, Manual, data sheet, scissors and two clamps 30 m 透析膜, 说明书, 数据表, 剪刀, 2只夹子	40224
		25 l 16	2.0	Dialysis membrane 5 m, Manual, data sheet 5 m 透析膜, 说明书, 数据表	40228-5
				Dialysis membrane 30 m, Manual, data sheet, scissors and two clamps 30 m 透析膜, 说明书, 数据表, 剪刀, 2只夹子	40228
		45 l 28.6	6.4	Dialysis membrane 5 m, Manual, data sheet 5 m 透析膜, 说明书, 数据表	40230-5
				Dialysis membrane 30 m, Manual, data sheet, scissors and two clamps 30 m 透析膜, 说明书, 数据表, 剪刀, 2只夹子	40230
T4	12-14 kDa	25 l 16	2.0	Dialysis membrane 5 m, Manual, data sheet 5 m 透析膜, 说明书, 数据表	40228-5
				Dialysis membrane 30 m, Manual, data sheet, scissors and two clamps 30 m 透析膜, 说明书, 数据表, 剪刀, 2只夹子	40228
		45 l 28.6	6.4	Dialysis membrane 5 m, Manual, data sheet 5 m 透析膜, 说明书, 数据表	40230-5
				Dialysis membrane 30 m, Manual, data sheet, scissors and two clamps 30 m 透析膜, 说明书, 数据表, 剪刀, 2只夹子	40230

Dialysis Accessories / 透析配件



Description 说明	Article number 产品型号
96-deep well plate (volume per well 2.2 ml) 96-孔板 (每孔 2.2ML)	40785
48-deep well plate (volume per well 5.0 ml) 48-孔板 (每孔 5.0ml)	40879
Float for XPRESS Micro Dialyzer (MD100/MD300) XPRESS 微透析器 (MD100/MD300) 浮动框	40750
Floater Small Round (9 qty MD100/MD300/ED300, 4 qty MD1000) 小型圆形浮动框 (9 qty MD100/MD300/ED300, 4 qty MD1000)	40795
Floater Medium Round (7 qty MD1000) 中型圆形浮动框 (7 qty MD1000)	40797
Sealing sheets for XPRESS Micro Dialyzer (MD100/MD300) XPRESS 微透析器 (MD100/MD300) 密封板	40751
Sealing sheet for deep well plate (included in Dialyzer kits) 多孔板密封板 (透析器试剂盒附件)	40949
Sample tubes for MD300 with screw caps, 10 pcs. MD300样品收集管, 螺丝帽, 10只	40799-10
Sample tubes for MD1000, with screw caps, 6 pcs. MD1000样品收集管, 螺丝帽, 6只	40880
Syringe filter PA, Ø 25 mm, pore size 0.2 microns, 100 pieces PA 针筒过滤器, 直径25 mm, 孔径0.2微米, 100只	40591
Syringe filter PA, Ø 25 mm, pore size 0.45 microns, 100 pieces CA 针筒过滤器, 直径25 mm, 孔径0.45 微米, 100只	40597
Syringe filter CA, Ø 30 mm, pore size 0.45 microns, 100 pieces CA 针筒过滤器, 直径30 mm, 孔径0.45 微米, 100只	40551
Syringe filter PTFE, Ø 15 mm, pore size 0.45 microns, 100 pieces PTFE 针筒过滤器, 直径15 mm, 孔径0.45 微米, 100只	40584
Syringe filter PTFE, Ø 25 mm, pore size 0.45 microns, 100 pieces PTFE 针筒过滤器, 直径25 mm, 孔径0.45 微米, 100只	40585

Cryo Labels / 冷冻标签



Dimensions [mm] 大小 (mm)	Sheets 数量 (张)	Article number 产品型号
24 x 20	20	80030
	50	80031
28 x 12	20	80020
	50	80021
33 x 14	20	80060
	50	80061
36 x 14	20	80010
	50	80011
67 x 25	20	80050
	50	80051
70 x 70	20	80040
	50	80041
size mixture 混合包装	18	80070

Labelling Software Barcode Forge scienova

Article 产品	Description 说明	Article number 产品型号
Barcode Forge scienova	workstation license for labeling software Barcode Forge scienova on USB-Stick including quick start guide USB-存储器, 含Barcode Forge scienova 标记软件的工作站许可证	81100
Startpackage Barcode Forge scienova	cryo labels 24 x 20, 20 sheets including software and quick start guide on usb-stick 低温标签, 24 x 20, 20张, 含软件和简易启动指南的存储器	80032
	cryo labels 28 x 12, 20 sheets including software and quick start guide on usb-stick 低温标签, 28 x 12, 20张, 含软件和简易启动指南的存储器	80022
	cryo labels 33 x 14, 20 sheets including software and quick start guide on usb-stick 低温标签, 33 x 14, 20张, 含软件和简易启动指南的存储器	80062
	cryo labels 36 x 14, 20 sheets including software and quick start guide on usb-stick 低温标签, 36 x 14, 20张, 含软件和简易启动指南的存储器	80012
	cryo labels 67 x 25, 20 sheets including software and quick start guide on usb-stick 低温标签, 67 x 25, 20张, 含软件和简易启动指南的存储器	80052
	cryo labels 70 x 70, 20 sheets including software and quick start guide on usb-stick 低温标签, 70 x 70, 20张, 含软件和简易启动指南的存储器	80042
	set consisting of 3 sheets of each dimension (total 18 sheets) including USB-Stick with Barcode Forge scienova and quick start guide	80072
Mixed package Barcode Forge scienova		

Barcode Forge
scienova 混合包装

一套, 每尺寸3张 (共18张), 含带scienova条形码和快速启动指南的存储器



请从www.scienova.com网
站下载免费试用版, 安
装, 试用



scienova GmbH • Spitzweidenweg 30 • 07743 Jena • Germany
Tel. +49 3641 504 586 • Fax +49 3641 504 587 • Info@scienova.com • www.scienova.com

© 2019 scienova GmbH | 2019 / 02